





## BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Sagittal Saw Systems

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®,  
Micro 100™ & Micro E®, All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3100C</b>	mm	41.5	14	0.40	0.61
	in	1.63	0.55	0.016	0.024
Replaces 5023-133, 5023-233 Replaces 5053-033, 5053-233					
					Coarse Straight Edge
<b>KM-3101</b>	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces 5023-138, 5023-238 Replaces 5053-038, 5053-238					
					Fine Straight Edge
<b>KM-3101C</b>	mm	25.7	9.5	0.40	0.61
	in	1.01	0.37	0.016	0.024
Replaces 5023-143, 5023-243 Replaces 5053-043, 5053-243					
					Coarse Straight Edge
<b>KM-3102</b>	mm	25.7	4.5	0.40	0.60
	in	1.01	0.18	0.016	0.024
Replaces 5023-139, 5023-239 Replaces 5053-039, 5053-239					
					Fine Straight Edge



CONMED® LINVATEC®  
Replacements





## BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Sagittal Saw Systems

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®,  
Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems







CONMED® LINVATEC®  
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3104</b>	mm	16.1	9.5	0.40	0.63
	in	0.63	0.37	0.016	0.025
Replaces N/A					
					Fine Straight Edge
<b>KM-3105</b>	mm	16.1	4.5	0.40	0.63
	in	0.63	0.18	0.016	0.025
Replaces N/A					
					Fine Straight Edge
<b>KM-3106</b>	mm	25.7	5.5	0.40	0.61
	in	1.01	0.22	0.016	0.024
Replaces 5023-132, 5023-232 Replaces 5053-032, 5053-232					
					Coarse Straight Edge
<b>KM-3107</b>	mm	34.3	6	0.40	0.60
	in	1.35	0.24	0.016	0.024
Replaces 5023-141, 5023-241 Replaces 5053-041, 5053-241					
					Fine Straight Edge

## BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3183</b>	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces 5023-183					
					Fine Straight Edge ACL Blade
<b>KM-3231</b>	mm	41.5	19.5	0.40	0.61
	in	1.63	0.77	0.016	0.024
Replaces 5023-131, 5023-231 Replaces 5053-031, 5053-231					
					Coarse Straight Edge
<b>KM-3237</b>	mm	25.7	14	0.40	0.61
	in	1.01	0.55	0.016	0.024
Replaces 5023-137, 5023-237 Replaces 5053-037, 5053-237					
					Coarse Straight Edge
<b>KM-3240</b>	mm	41.5	9.5	0.40	0.60
	in	1.63	0.37	0.016	0.024
Replaces 5023-140, 5023-240 Replaces 5053-040, 5053-240					
					Fine Straight Edge Tapered

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems





**CONMED® LINVATEC®**  
Replacements

## BUSA® Small Bone Sagittal Saw Blades - Replacements - Linvatec®





Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3179SN</b>	mm	42.5	54	0.40	0.69
	in	1.67	2.13	0.016	0.027
Replaces N/A					
 <p style="text-align: right;">Coarse Opposing Teeth Non-Sterile</p>					
<b>KM-3231SN</b>	mm	41.9	19	0.40	0.71
	in	1.65	0.75	0.016	0.028
Replaces N/A					
 <p style="text-align: right;">Non-Sterile</p>					

## BUSA® Small Bone Sagittal Saw Blades Graft Harvesting - Replacements - Linvatec®

Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;  
All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3405</b>	mm	20	16.5	0.40	5.0
	in	0.79	0.65	0.016	0.197
Replaces 5023-405					
 <p style="text-align: right;">Graft Harvesting Blade</p>					
<b>KM-3406</b>	mm	20	16.5	0.40	6.0
	in	0.79	0.65	0.016	0.236
Replaces 5023-406					
 <p style="text-align: right;">Graft Harvesting Blade</p>					
<b>KM-3407</b>	mm	20	16.5	0.40	7.0
	in	0.79	0.65	0.016	0.276
Replaces 5023-407					
 <p style="text-align: right;">Graft Harvesting Blade</p>					
<b>KM-3408</b>	mm	20	16.5	0.40	8.0
	in	0.79	0.65	0.016	0.315
Replaces 5023-408					
 <p style="text-align: right;">Graft Harvesting Blade</p>					



Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®, All BUSA® Small Bone Sagittal Saw Systems



**CONMED® LINVATEC®**  
Replacements

## BUSA® Small Bone Sagittal Saw Blades Graft Harvesting - Replacements - Linvatec®





Compatible with Linvatec® Sagittal Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Sagittal Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3409</b>	mm	20	16.5	0.40	9.0
	in	0.79	0.65	0.016	0.354
Replaces 5023-409					
 <p style="text-align: right;">Graft Harvesting Blade</p>					
<b>KM-3410</b>	mm	20	16.5	0.40	10.0
	in	0.79	0.65	0.016	0.394
Replaces 5023-410					
 <p style="text-align: right;">Graft Harvesting Blade</p>					

## BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Oscillating Saw Systems

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®,  
Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3000</b>	mm	18.5	5.5	0.40	0.60
	in	0.73	0.22	0.016	0.024
Replaces 5023-162, 5023-262 Replaces 5053-062, 5053-262					
 <div style="text-align: right;">Fine Straight Edge</div>					
<b>KM-3001</b>	mm	28.5	5.5	0.40	0.60
	in	1.12	0.22	0.016	0.024
Replaces 5023-164, 5023-264 Replaces 5053-064, 5053-264					
 <div style="text-align: right;">Fine Straight Edge</div>					
<b>KM-3002</b>	mm	34.7	10	0.40	0.60
	in	1.37	0.39	0.016	0.024
Replaces 5023-169, 5023-269 Replaces 5053-073, 5053-273					
 <div style="text-align: right;">Fine Straight Edge</div>					
<b>KM-3003</b>	mm	22	5.5	0.40	0.60
	in	0.87	0.22	0.016	0.024
Replaces N/A					
 <div style="text-align: right;">Fine Straight Edge</div>					







CONMED® LINVATEC®  
Replacements

## BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Oscillating Saw Systems

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®,  
Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems



CONMED® LINVATEC®  
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3004</b>	mm	27.9	9.5	0.40	0.60
	in	1.10	0.37	0.016	0.024
Replaces N/A					
					Fine Straight Edge
<b>KM-3005</b>	mm	18.5	10	0.40	0.60
	in	0.73	0.39	0.016	0.024
Replaces 5023-163, 5023-263 Replaces 5053-063, 5053-263					
					Fine Straight Edge
<b>KM-3006</b>	mm	12	5.5	0.40	0.60
	in	0.47	0.22	0.016	0.024
Replaces 5023-160, 5023-260 Replaces 5053-060, 5053-260					
					Fine Straight Edge
<b>KM-3008</b>	mm	28.5	18.5	0.40	0.60
	in	1.12	0.73	0.016	0.024
Replaces 5023-166, 5023-266 Replaces 5053-066, 5053-266					
					Fine Radial Edge



## BUSA® Small Bone Oscillating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®;  
All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3261</b>	mm	12	10.0	0.40	0.60
	in	0.47	0.39	0.016	0.024
Replaces 5023-161, 5023-261 Replaces 5053-061, 5053-261					
 <div style="text-align: right;">Fine Straight Edge</div>					
<b>KM-3265</b>	mm	28.5	10.00	0.40	0.60
	in	1.12	0.39	0.016	0.024
Replaces 5023-165, 5023-265 Replaces 5053-065, 5053-265					
 <div style="text-align: right;">Fine Radial Edge</div>					





Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro I 00™ & Micro E®, All BUSA® Small Bone Oscillating Saw Systems



**CONMED® LINVATEC®**  
Replacements




## BUSA® Small Bone Oscillating Saw Blades Crescentic - Replacements - Linvatec®

Compatible with Linvatec® Oscillating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Oscillating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3170</b>	mm	12	9	0.40	0.60
	in	0.47	0.36	0.016	0.024
Replaces 5023-170, 5023-270 Replaces 5053-170, 5053-270					
 <div style="text-align: right;"><i>Fine Crescentic Blade</i></div>					
<b>KM-3171</b>	mm	25.5	13.8	0.40	0.60
	in	1.00	0.55	0.016	0.024
Replaces 5023-171, 5023-271 Replaces 5053-171, 5053-271					
 <div style="text-align: right;"><i>Fine Crescentic Blade</i></div>					
<b>KM-3172</b>	mm	32	18	0.40	0.60
	in	1.26	0.71	0.016	0.024
Replaces 5023-172, 5023-272 Replaces 5053-172, 5053-272					
 <div style="text-align: right;"><i>Fine Crescentic Blade</i></div>					
<b>KM-3176</b>	mm	32	13.8	0.40	0.60
	in	1.26	0.55	0.016	0.024
Replaces 5023-176, 5023-276 Replaces 5053-176, 5053-276					
 <div style="text-align: right;"><i>Fine Crescentic Blade</i></div>					

## BUSA® Small Bone Reciprocating Saw Blades - Replacements - Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®; All BUSA® Small Bone Reciprocating Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-060</b>	mm	25.5	7.6	0.40	0.64
	in	1.00	0.30	0.016	0.025
Replaces 5053-020, 5053-220					
					Flat Shank
<b>KM-070</b>	mm	38	7.6	0.40	0.64
	in	1.50	0.30	0.016	0.025
Replaces 5053-021, 5053-221					
					Flat Shank
<b>KM-54</b>	mm	23.2	7.4	0.40	0.60
	in	0.913	0.29	0.016	0.024
Replaces 5053-054, 5053-254					
					Round Shank

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®, All BUSA® Small Bone Reciprocating Saw Systems





CONMED® LINVATEC®  
Replacements

## BUSA® Small Bone Reciprocating Rasps - Replacements - Linvatec®

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®, Micro 100™ & Micro E®;  
All BUSA® Small Bone Reciprocating Saw Systems

Compatible with Linvatec® Reciprocating Saw Systems: MicroPower®, MicroChoice®,  
Micro 100™ & Micro E®; All BUSA® Small Bone Reciprocating Saw Systems

CONMED® LINVATEC®  
Replacements

Order Number	Unit	Length	Max. Width
<b>KM-80</b>	mm	11	5.5
	in	0.43	0.22
Replaces 5053-080			
			Small Tear Round Shank
<b>KM-81</b>	mm	13	6.9
	in	0.51	0.27
Replaces 5053-081			
			Large Tear Round Shank
<b>KM-82</b>	mm	11	5.5
	in	0.43	0.22
Replaces 5053-082			
			Small Tear Cross-Cut Round Shank
<b>KM-83</b>	mm	13	6.9
	in	0.51	0.27
Replaces 5053-083			
			Large Tear Cross-Cut Round Shank