







BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0901	Oval - Carbide	mm	44.5	4.0	7.9	8
		in	1.75	0.16	0.31	
Replaces 277-10-62						
						1 per Pack Sterile
14-C0905	Oval - Carbide	mm	44.5	5.5	10.3	10
		in	1.75	0.22	0.41	
Replaces N/A						
						1 per Pack Sterile
16-C0901	Oval - Carbide	mm	54.5	4.0	7.9	8
		in	2.15	0.16	0.31	
Replaces N/A						
						1 per Pack Sterile
14-C0103	Round - Carbide	mm	44.5	1.5	1.5	6
		in	1.75	0.06	0.06	
Replaces N/A						
						1 per Pack Sterile
14-C0104	Round - Carbide	mm	44.5	1.8	1.8	6
		in	1.75	0.07	0.07	
Replaces 277-10-418						
						1 per Pack Sterile
14-C0105	Round - Carbide	mm	44.5	2.0	2.0	6
		in	1.75	0.08	0.08	
Replaces 277-10-225						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0106	Round - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces 277-10-423						
						8 Shannon 1 per Pack Sterile
14-C0107	Round - Carbide	mm	44.5	3.0	3.0	6
		in	1.75	0.12	0.12	
Replaces 277-10-230						
						1 per Pack Sterile
14-C0108	Round - Carbide	mm	44.5	4.0	4.0	8
		in	1.75	0.16	0.16	
Replaces 277-10-240						
						1 per Pack Sterile
14-C0109	Round - Carbide	mm	44.5	5.0	5.0	8
		in	1.75	0.20	0.20	
Replaces 277-10-250						
						1 per Pack Sterile
14-C0110	Round - Carbide	mm	44.5	6.0	6.0	8
		in	1.75	0.24	0.24	
Replaces N/A						
						1 per Pack Sterile
16-C0105	Round - Carbide	mm	54.5	2.0	2.0	6
		in	2.15	0.08	0.08	
Replaces N/A						
						1 per Pack Sterile

Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-C0107	Round - Carbide	mm	54.5	3.0	3.0	6
		in	2.15	0.12	0.12	
Replaces N/A						
						1 per Pack Sterile
16-C0108	Round - Carbide	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
Replaces N/A						
						1 per Pack Sterile
14-C0308	Straight - Cross-Cut - Carbide	mm	44.5	1.6	4.5	6
		in	1.75	0.06	0.17	
Replaces 1607-2-103						
						560 1 per Pack Sterile
14-C0405	Straight - Cross-Cut - Carbide	mm	44.5	1.6	5.9	6
		in	1.75	0.06	0.23	
Replaces N/A						
						1 per Pack Sterile
14-C0606	Tapered - Cross-Cut - Carbide	mm	44.5	1.2	4.1	6
		in	1.75	0.05	0.16	
Replaces 277-10-212						
						1 per Pack Sterile
14-C0608	Tapered - Cross-Cut - Carbide	mm	44.5	1.6	4.3	6
		in	1.75	0.06	0.17	
Replaces 277-10-216						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0616	Tapered - Cross-Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces 277-10-221						
						1 per Pack Sterile
16-D0129	Round - Diamond	mm	54.5	1.5	1.5	Fine
		in	2.15	0.06	0.06	
Replaces 1608-6-85						
						1 per Pack Sterile
16-D0130	Round - Diamond	mm	54.5	2.0	2.0	Fine
		in	2.15	0.08	0.08	
Replaces 1608-6-87						
						1 per Pack Sterile
16-D0131	Round - Diamond	mm	54.5	2.4	2.4	Fine
		in	2.15	0.09	0.09	
Replaces 1608-6-89						
						1 per Pack Sterile
16-D0132	Round - Diamond	mm	54.5	3.0	3.0	Medium
		in	2.15	0.11	0.12	
Replaces 1608-6-91						
						1 per Pack Sterile
16-D0134	Round - Diamond	mm	54.5	3.5	3.5	Medium
		in	2.15	0.14	0.14	
Replaces 1608-6-93						
						1 per Pack Sterile

Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-D0135	Round - Diamond	mm	54.5	4.0	4.0	Medium
		in	2.15	0.16	0.16	
Replaces 1608-6-95						
						1 per Pack Sterile
16-D0137	Round - Diamond	mm	54.5	5.0	5.0	Medium
		in	2.15	0.20	0.20	
Replaces 1608-6-99						
						1 per Pack Sterile
14-S0402	Cylinder - Stainless Steel	mm	44.5	6.1	11.1	14
		in	1.75	0.24	0.44	
Replaces 1607-2-1						
						1 per Pack Sterile
16-S0402	Cylinder - Stainless Steel	mm	54.5	6.1	11.1	14
		in	2.15	0.24	0.44	
Replaces 1608-2-1						
						1 per Pack Sterile
15-S1101	Flame - Stainless Steel	mm	48.0	6.0	12.7	14
		in	1.89	0.24	0.50	
Replaces 1607-2-47						
						1 per Pack Sterile
14-S0201	Inverted Cone - Stainless Steel	mm	44.5	3.1	3.3	10
		in	1.75	0.12	0.13	
Replaces 1607-2-57						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-S1804	Lindemann - Stainless Steel	mm	44.5	1.6	9.0	4
		in	1.75	0.06	0.35	
<i>Replaces N/A</i>						
						<i>1 per Pack Sterile</i>
16-S1805	Lindemann - Stainless Steel	mm	54.5	1.6	9.0	4
		in	2.15	0.06	0.35	
<i>Replaces N/A</i>						
						<i>1 per Pack Sterile</i>
16-S1806	Lindemann - Stainless Steel	mm	54.5	2.1	10.0	4
		in	2.15	0.08	0.39	
<i>Replaces N/A</i>						
						<i>1 per Pack Sterile</i>
16-S1801	Lindemann - Stainless Steel	mm	54.5	2.3	9.5	4
		in	2.15	0.09	0.37	
<i>Replaces 1608-2-21</i>						
						<i>1 per Pack Sterile</i>
19-S1802	Lindemann - Stainless Steel	mm	66.7	2.2	24.3	4
		in	2.63	0.09	0.96	
<i>Replaces 1608-2-43</i>						
						<i>1 per Pack Sterile</i>
23-S1809	Lindemann - Stainless Steel	mm	70.0	2.3	22.0	4
		in	2.76	0.09	0.87	
<i>Replaces N/A</i>						
						<i>1 per Pack Sterile</i>

Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
22-S1808	Lindemann - Stainless Steel	mm	75.0	2.3	35.0	4
		in	2.95	0.09	1.38	
Replaces N/A						
						1 per Pack Sterile
14-S0903	Oval - Stainless Steel	mm	44.5	4.0	9.5	10
		in	1.75	0.16	0.37	
Replaces 1607-2-35						
						1 per Pack Sterile
14-S0909	Oval - Stainless Steel	mm	44.5	6.0	9.5	14
		in	1.75	0.24	0.37	
Replaces 1607-2-3						
						1 per Pack Sterile
16-S0903	Oval - Stainless Steel	mm	54.5	4.0	9.5	10
		in	2.15	0.16	0.37	
Replaces 1608-2-35						
						1 per Pack Sterile
16-S0909	Oval - Stainless Steel	mm	54.5	6.0	9.5	14
		in	2.15	0.24	0.37	
Replaces 1608-2-3						
						1 per Pack Sterile
14-S1201	Pear - Stainless Steel	mm	44.5	6.1	10.0	14
		in	1.75	0.24	0.39	
Replaces 1607-2-5						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S1201	Pear - Stainless Steel	mm	54.5	6.1	10.0	14
		in	2.15	0.24	0.39	
Replaces 1608-2-5						
						1 per Pack Sterile
14-S0117	Round - Stainless Steel	mm	44.5	3.1	3.1	8
		in	1.75	0.12	0.12	
Replaces N/A						
						1 per Pack Sterile
14-S0124	Round - Stainless Steel	mm	44.5	6.0	6.0	16
		in	1.75	0.24	0.24	
Replaces 1607-2-15						
						1 per Pack Sterile
14-S0125	Round - Stainless Steel	mm	44.5	7.0	7.0	16
		in	1.75	0.28	0.28	
Replaces 1607-2-17						
						1 per Pack Sterile
14-S0127	Round - Stainless Steel	mm	44.5	8.0	8.0	16
		in	1.75	0.31	0.31	
Replaces 1607-2-41						
						1 per Pack Sterile
16-S0112	Round - Stainless Steel	mm	54.5	1.4	1.4	16
		in	2.15	0.055	0.055	
Replaces 1608-6-157						
						1 per Pack Sterile

Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0103	Round - Stainless Steel	mm	54.5	1.5	1.5	6
		in	2.15	0.06	0.06	
Replaces 1608-2-7						
						1 per Pack Sterile
16-S0106	Round - Stainless Steel	mm	54.5	2.3	2.3	6
		in	2.15	0.09	0.09	
Replaces 1608-6-147						
						1 per Pack Sterile
16-S0115	Round - Stainless Steel	mm	54.5	2.4	2.4	8
		in	2.15	0.095	0.095	
Replaces 1608-2-9						
						1 per Pack Sterile
16-S0107	Round - Stainless Steel	mm	54.5	3.2	3.2	6
		in	2.15	0.126	0.126	
Replaces 1608-6-149						
						1 per Pack Sterile
16-S0117	Round - Stainless Steel	mm	54.5	3.2	3.2	8
		in	2.15	0.125	0.125	
Replaces 1608-2-11						
						1 per Pack Sterile
16-S0116	Round - Stainless Steel	mm	54.5	3.0	3.0	20
		in	2.15	0.12	0.12	
Replaces 1608-6-163						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0119	Round - Stainless Steel	mm	54.5	4.0	4.0	6
		in	2.15	0.16	0.16	
<i>Replaces 1608-6-151</i>						
						<i>1 per Pack Sterile</i>
16-S0108	Round - Stainless Steel	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
<i>Replaces 1608-6-137</i>						
						<i>1 per Pack Sterile</i>
16-S0120	Round - Stainless Steel	mm	54.5	4.0	4.0	12
		in	2.15	0.16	0.16	
<i>Replaces 1608-2-55</i>						
						<i>1 per Pack Sterile</i>
16-S0118	Round - Stainless Steel	mm	54.5	4.0	4.0	24
		in	2.15	0.16	0.16	
<i>Replaces 1608-6-165</i>						
						<i>1 per Pack Sterile</i>
16-S0121	Round - Stainless Steel	mm	54.5	4.6	4.6	6
		in	2.15	0.18	0.18	
<i>Replaces 1608-6-153</i>						
						<i>1 per Pack Sterile</i>
16-S0122	Round - Stainless Steel	mm	54.5	4.6	4.6	8
		in	2.15	0.18	0.18	
<i>Replaces 1608-6-139</i>						
						<i>1 per Pack Sterile</i>

Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0123	Round - Stainless Steel	mm	54.5	4.5	4.5	12
		in	2.15	0.177	0.177	
Replaces 1608-2-13						
						1 per Pack Sterile
16-S0109	Round - Stainless Steel	mm	54.5	5.1	5.1	8
		in	2.15	0.20	0.20	
Replaces 1608-6-141						
						1 per Pack Sterile
16-S0110	Round - Stainless Steel	mm	54.5	6.1	6.1	8
		in	2.15	0.24	0.24	
Replaces 1608-6-143						
						1 per Pack Sterile
16-S0124	Round - Stainless Steel	mm	54.5	6.0	6.0	16
		in	2.15	0.24	0.24	
Replaces 1608-2-15						
						1 per Pack Sterile
16-S0126	Round - Stainless Steel	mm	54.5	7.1	7.1	8
		in	2.15	0.28	0.28	
Replaces 1608-6-145						
						1 per Pack Sterile
16-S0125	Round - Stainless Steel	mm	54.5	7.0	7.0	16
		in	2.15	0.28	0.28	
Replaces 1608-2-17						
						1 per Pack Sterile

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S1401	Shannon Style - Tapered - Stainless Steel	mm	54.5	2.1	11.1	3
		in	2.15	0.08	0.44	
<i>Replaces 1608-2-101</i>						
						<i>1 per Pack Sterile</i>
16-S0313	Straight - Stainless Steel	mm	54.5	1.5	11.7	3
		in	2.15	0.06	0.46	
<i>Replaces 1608-2-51</i>						
						<i>1 per Pack Sterile</i>
16-S0621	Tapered - Cross-Cut - Stainless Steel	mm	54.5	3.1	11.1	8
		in	2.15	0.12	0.44	
<i>Replaces 1608-6-167</i>						
						<i>1 per Pack Sterile</i>
16-S1601	Swanson Reamer - Stainless Steel	mm	56.5	2.0	10.5	5
		in	2.23	0.08	0.41	
<i>Replaces 1608-2-131</i>						
						<i>1 per Pack Sterile</i>
17-S1602	Swanson Reamer - Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
<i>Replaces 1608-2-133</i>						
						<i>1 per Pack Sterile</i>
17-S1603	Swanson Reamer - Stainless Steel	mm	62.0	4.0	13.5	6
		in	2.44	0.16	0.53	
<i>Replaces 1608-2-135</i>						
						<i>1 per Pack Sterile</i>





Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

STRYKER®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery J-Notch - Replacements - Stryker®

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™ & Command 2







Compatible with Stryker® Micro Drill Systems:
RemB™, CORE, TPS™ & Command 2

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S1702	Wire Pass Drill - Stainless Steel	mm	54.5	1.2	10.0	2
		in	2.15	0.047	0.39	
<i>Replaces 1608-2-71</i>						
						<i>1 per Pack Sterile</i>
18-S1701	Wire Pass Drill - Stainless Steel	mm	64.5	1.0	10.0	2
		in	2.54	0.039	0.39	
<i>Replaces 1608-2-63</i>						
						<i>1 per Pack Sterile</i>
18-S1703	Wire Pass Drill - Stainless Steel	mm	65.5	1.5	6.9	2
		in	2.57	0.059	0.27	
<i>Replaces 1608-2-59</i>						
						<i>1 per Pack Sterile</i>
19-S1705	Wire Pass Drill - Stainless Steel	mm	66.7	2.0	6.4	2
		in	2.63	0.079	0.25	
<i>Replaces 1608-2-61</i>						
						<i>1 per Pack Sterile</i>

STRYKER®
Replacements

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KM946N	Round - Carbide	mm	44.5	1.8	1.8	6
		in	1.75	0.07	0.07	
Replaces N/A						
						5 per Pack Non-Sterile
KMS946N	Round - Carbide	mm	44.5	1.8	1.8	6
		in	1.75	0.07	0.07	
Replaces 2296-100-418						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM948N	Round - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces N/A						
						5 per Pack Non-Sterile
KMS948N	Round - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces 2296-100-423						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM956N	Round - Carbide	mm	58.4	1.8	1.8	6
		in	2.30	0.07	0.07	
Replaces N/A						
						5 per Pack Non-Sterile
KMS956N	Round - Carbide	mm	58.4	1.8	1.8	6
		in	2.30	0.07	0.07	
Replaces 2296-101-418						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)

Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max







STRYKER®
Replacements

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max







Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KM958N	Round - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces N/A						
						5 per Pack Non-Sterile
KMS958N	Round - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces 2296-101-423						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM948SN	Round - Super Sharp - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces N/A						
						5 per Pack Non-Sterile
KM958SN	Round - Super Sharp - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces N/A						
						5 per Pack Non-Sterile
KMS958SN	Round - Super Sharp - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces N/A						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM94701N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	1.2	4.2	6
		in	1.75	0.05	0.17	
Replaces N/A						
						5 per Pack Non-Sterile

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS94701N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	1.2	4.2	6
		in	1.75	0.05	0.17	
Replaces 2296-100-212						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM94702N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	1.6	4.4	6
		in	1.75	0.06	0.17	
Replaces N/A						
						5 per Pack Non-Sterile
KMS94702N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	1.6	4.4	6
		in	1.75	0.06	0.17	
Replaces 2296-100-216						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM94703N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces N/A						
						5 per Pack Non-Sterile
KMS94703N	Tapered Fissure - Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces 2296-100-221						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM95701N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	1.2	4.2	6
		in	2.30	0.05	0.17	
Replaces N/A						
						5 per Pack Non-Sterile







Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max

STRYKER®
Replacements

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max







Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS95701N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	1.2	4.2	6
		in	2.30	0.05	0.17	
Replaces 2296-101-212						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM95702N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	1.6	4.4	6
		in	2.30	0.06	0.17	
Replaces N/A						
						5 per Pack Non-Sterile
KMS95702N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	1.6	4.4	6
		in	2.30	0.06	0.17	
Replaces 2296-101-216						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM95703N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	2.1	4.8	6
		in	2.30	0.08	0.19	
Replaces N/A						
						5 per Pack Non-Sterile
KMS95703N	Tapered Fissure - Cross Cut - Carbide	mm	58.4	2.1	4.8	6
		in	2.30	0.08	0.19	
Replaces 2296-101-221						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM94702LN	Tapered Fissure - Long Cross Cut - Carbide	mm	44.5	1.6	6.0	6
		in	1.75	0.06	0.24	
Replaces N/A						
						5 per Pack Non-Sterile

STRYKER®
Replacements

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS94702LN	Tapered Fissure - Long Cross Cut - Carbide	mm	44.5	1.6	6.0	6
		in	1.75	0.06	0.24	
Replaces 2296-100-316						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM94703LN	Tapered Fissure - Long Cross Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces N/A						
						5 per Pack Non-Sterile
KMS94703LN	Tapered Fissure - Long Cross Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces 2296-100-321						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM95702LN	Tapered Fissure - Long Cross Cut - Carbide	mm	58.4	1.6	6.0	6
		in	2.30	0.06	0.24	
Replaces N/A						
						5 per Pack Non-Sterile
KMS95702LN	Tapered Fissure - Long Cross Cut - Carbide	mm	58.4	1.6	6.0	6
		in	2.30	0.06	0.24	
Replaces 2296-101-316						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM95703LN	Tapered Fissure - Long Cross Cut - Carbide	mm	58.4	2.1	7.5	6
		in	2.30	0.08	0.30	
Replaces N/A						
						5 per Pack Non-Sterile

Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max







STRYKER®
Replacements

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max







Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max

STRYKER®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS95703LN	Tapered Fissure - Long Cross Cut - Carbide	mm	58.4	2.1	7.5	6
		in	2.30	0.08	0.30	
Replaces 2296-101-321						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM941703LN	Tapered Fissure - Long Round End - Cross-Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces N/A						
						5 per Pack Non-Sterile
KMS941703LN	Tapered Fissure - Long Round End - Cross-Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces 2296-100-421						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM951703LN	Tapered Fissure - Long Round End - Cross-Cut - Carbide	mm	58.4	2.1	7.5	6
		in	2.30	0.08	0.30	
Replaces N/A						
						5 per Pack Non-Sterile
KMS951703LN	Tapered Fissure - Long Round End - Cross-Cut - Carbide	mm	58.4	2.1	7.5	6
		in	2.30	0.08	0.30	
Replaces 2296-101-421						
						5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)
KM941703SN	Tapered Super Sharp Fissure - Long Round End - Cross-Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces N/A						
						5 per Pack Non-Sterile

BUSA® Micro Burs - Oral Surgery Flat Shank - Replacements - Stryker®

Compatible with Stryker® Impaction Drills: 100K, CMD, CORE & Oral Max

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS941703SN	Tapered Super Sharp Fissure - Long Round End - Cross-Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces N/A						
 <p style="text-align: right;">5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)</p>						
KM951703SN	Tapered Super Sharp Fissure - Long Round End - Cross-Cut - Carbide	mm	58.4	2.1	4.8	6
		in	2.30	0.08	0.19	
Replaces N/A						
 <p style="text-align: right;">5 per Pack Non-Sterile</p>						
KMS951703SN	Tapered Super Sharp Fissure - Long Round End - Cross-Cut - Carbide	mm	58.4	2.1	4.8	6
		in	2.30	0.08	0.19	
Replaces N/A						
 <p style="text-align: right;">5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)</p>						
KM94901N	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces 2296-100-026						
 <p style="text-align: right;">5 per Pack Non-Sterile</p>						
KM95901N	Oval - Carbide	mm	58.4	4.0	8.0	8
		in	2.30	0.16	0.31	
Replaces N/A						
 <p style="text-align: right;">5 per Pack Non-Sterile</p>						
KMS95901N	Oval - Carbide	mm	58.4	4.0	8.0	8
		in	2.30	0.16	0.31	
Replaces 2296-101-026						
 <p style="text-align: right;">5 per Pack Sterile (Individual Sterile Packs with Lot Numbers)</p>						

Compatible with Stryker® Impaction Drills:
100K, CMD, CORE & Oral Max

STRYKER®
Replacements