

# SHUKLA™ Medical

Innovators of  
**Xtract-All®**

## Universal Stripped and Seized Screw Extraction System


### Removes:

- M2 to M8 Screws
- Cold Weld (fused) Screws  
(Capable of applying over 250 in-lbs of Torque)

### Attachments include:

- Extractors
- 3 Extensions
- Breaker Bar
- Ratcheting T Handle

Have the instrument you need  
when you need it.




Contact us at: **888-4shukla (888-474-8552)**  
or **customer.service@ShuklaMedical.com**  
See our website at: **www.ShuklaMedical.com**



# Universal Stripped and Seized Screw Extraction System

**SHUKLA**<sup>TM</sup>  
Medical



(MKS 1019)

## STRIPPED AND SEIZED BONE SCREW EXTRACTION SYSTEM

No.	Part Number	Description	Units	Single Use
1	TP347-M01	Extractor tip for 2.0-2.5 mm Screw Head Recess	4	✓
2	TP348-M01	Extractor tip for 2.5-3.5 mm Screw Head Recess	4	✓
3	TP349-M01	Extractor tip for 3.5-4.5 mm Screw Head Recess	4	✓
4	TP350-M01	Extractor tip for 4.5-5.5 mm Screw Head Recess	4	✓
5	TP351-M01	Extractor tip for 5.5-6.5 mm Screw Head Recess	4	✓
6	TP352-M01	Extractor tip for 6.0-7.0 mm Screw Head Recess	4	✓
7	MMI1152-M01	4" Extension Driver	1	
8	MMI1153-M01	6" Extension Driver	1	
9	MMI1154-M01	8" Extension Driver	1	
10	CA137	Case, Stripped and Seized Screw Removal System	1	
11	CA138	Lid, Stripped and Seized Screw Removal System	1	
12	CA139	Caddy, Stripped and Seized Screw Removal System	1	
13	CA140	Caddy lid, Stripped and Seized Screw Removal System	1	
14	HD233-M01	Asymmetric Cannulated Ratcheting T-Handle	1	
15	HD239-M01	Breaker Bar	1	

