

SHUKLA™ Medical

Innovators of
Xtract-All®

Universal Broken Screw Extraction System

Removes:

- M2 to M8 Size Screws

Attachments include:

- Extractors
- 3 Extensions
- Jacobs Chuck Attachment
- Ratcheting T Handle




**Have the instrument you need
when you need it.**

Contact us at: **888-4shukla (888-474-8552)**
or **customer.service@ShuklaMedical.com**
See our website at: **www.ShuklaMedical.com**



Universal Broken Screw Extraction System

SHUKLATM
Medical



(MKS 1020)

BROKEN BONE SCREW EXTRACTION SYSTEM

No.	Part Number	Description	Units	Single Use
1	TP341-M01	Extractor tip for screw with 2.0-2.5 mm Thread	4	✓
2	TP342-M01	Extractor tip for screw with 2.5-4.0 mm Thread	4	✓
3	TP343-M01	Extractor tip for screw with 3.5-5.0 mm Thread	4	✓
4	TP344-M01	Extractor tip for screw with 4.5-6.0 mm Thread	4	✓
5	TP345-M01	Extractor tip for screw with 5.5-7.0 mm Thread	4	✓
6	TP346-M01	Extractor tip for screw with 7.0-8.0 mm Thread	4	✓
7	MMI1152-M01	4" Extension Driver	1	
8	MMI1153-M01	6" Extension Driver	1	
9	MMI1154-M01	8" Extension Driver	1	
10	CA141	Case, Broken Screw Extraction System	1	
11	CA142	Lid, Broken Screw Extraction System	1	
12	CA143	Caddy, Broken Screw Extraction System	1	
13	CA144	Caddy lid, Broken Screw Extraction System	1	
14	HD233-M01	Asymmetric Cannulated Ratcheting T-Handle	1	
15	AD120-M01	Power Drill Adapter	1	

