





BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012600	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
<i>No Tooth Offset</i>					
					
PMX-012601	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
<i>No Tooth Offset</i>					
					
PMX-012602	mm	90	12.5	0.99	0.99
	in	3.54	0.49	0.039	0.039
<i>No Tooth Offset</i>					
					
PMX-012603	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
<i>No Tooth Offset</i>					
					

Compatible with All BUSA®
Sagittal/Oscillator Saw Systems







BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems





Compatible with All BUSA®
Sagittal/Oscillator Saw Systems

BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012604	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
					No Tooth Offset
PMX-012610	mm	100	19	1.27	1.27
	in	3.94	0.75	0.050	0.050
					No Tooth Offset
PMX-012611	mm	75	19	0.89	0.89
	in	2.95	0.75	0.035	0.035
					No Tooth Offset
PMX-012612	mm	90	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
					No Tooth Offset

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012613	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
					No Tooth Offset
PMX-012614	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
					No Tooth Offset
PMX-012615	mm	100	19	1.19	1.19
	in	3.94	0.75	0.047	0.047
					No Tooth Offset
PMX-012616	mm	90	19	0.99	0.99
	in	3.54	0.75	0.039	0.039
					No Tooth Offset

Compatible with All BUSA®
Sagittal/Oscillator Saw Systems







BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems





Compatible with All BUSA®
Sagittal/Oscillator Saw Systems

BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012620	mm	100	25	1.27	1.27
	in	3.94	0.98	0.050	0.050
					No Tooth Offset
PMX-012621	mm	75	25	0.89	0.89
	in	2.95	0.98	0.035	0.035
					No Tooth Offset
PMX-012622	mm	90	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
					No Tooth Offset
PMX-012623	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
					No Tooth Offset

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012624	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
					No Tooth Offset
PMX-012625	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
					No Tooth Offset
PMX-012626	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
					No Tooth Offset
PMX-012627	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
					No Tooth Offset

Compatible with All BUSA®
Sagittal/Oscillator Saw Systems







BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems




Compatible with All BUSA®
Sagittal/Oscillator Saw Systems

BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012628	mm	95	19	1.37	1.37
	in	3.74	0.75	0.054	0.054
					No Tooth Offset
PMX-012629	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
					No Tooth Offset
					Flared Blade
PMX-012630	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
					No Tooth Offset
					Flared Blade
PMX-012631	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
					No Tooth Offset
					Flared Blade

BUSA® Large Bone Saw Blades - Sagittal/Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with All BUSA® Sagittal/Oscillator Saw Systems

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
PMX-012632	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
					
PMX-012633	mm	90	31	0.99	0.99
	in	3.54	1.22	0.039	0.039
					
PMX-012634	mm	90	31	1.37	1.37
	in	3.54	1.22	0.054	0.054
					

Compatible with All BUSA®
Sagittal/Oscillator Saw Systems



BUSA®
Proprietary