

# BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Compatible with DePuy® Synthes®  
Large Bone Oscillating Saw System







DEPUY® SYNTHES®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>K5-1270-50F</b>	mm	70	12.5	1.27	1.27
	in	2.76	0.49	0.050	0.050
<i>Replaces N/A</i>					
					<p>No Tooth Offset</p> <p>Flared Blade</p>
<b>K5-1290-47</b>	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
<i>Replaces 05.003.1045</i>					
					No Tooth Offset
<b>K5-1290-50</b>	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
<i>Replaces 05.003.1065</i>					
					No Tooth Offset
<b>K5-1290-54</b>	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
<i>Replaces 05.003.1075</i>					
					No Tooth Offset

# BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>K5-1990-39</b>	mm	90	19	1.00	1.00
	in	3.54	0.75	0.039	0.039
<i>Replaces 05.003.1125</i>					
					No Tooth Offset
<b>K5-1990-47</b>	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
<i>Replaces 05.003.1145</i>					
					No Tooth Offset
<b>K5-1990-50</b>	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
<i>Replaces 05.003.1165</i>					
					No Tooth Offset
<b>K5-1990-54</b>	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
<i>Replaces 05.003.1175</i>					
					No Tooth Offset

Compatible with DePuy® Synthes®  
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

# BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Compatible with DePuy® Synthes®  
Large Bone Oscillating Saw System







DEPUY® SYNTHES®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>K5-2190-39F</b>	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
<i>Replaces 05.002.306S</i>					
					No Tooth Offset
Flared Blade					
<b>K5-2190-47F</b>	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
<i>Replaces 05.002.302S</i>					
					No Tooth Offset
Flared Blade					
<b>K5-2190-50F</b>	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
<i>Replaces 05.002.303S</i>					
					No Tooth Offset
Flared Blade					
<b>K5-2190-54F</b>	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
<i>Replaces 05.002.304S</i>					
					No Tooth Offset
Flared Blade					

# BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>K5-2590-39</b>	mm	90	25	1.00	1.00
	in	3.54	0.98	0.039	0.039
<i>Replaces 05.003.132S</i>					
					<p>No Tooth Offset</p> <p>Flared Blade</p>
<b>K5-2590-47</b>	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
<i>Replaces 05.003.134S</i>					
					No Tooth Offset
<b>K5-2590-50</b>	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
<i>Replaces 05.003.136S</i>					
					No Tooth Offset
<b>K5-2590-54</b>	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
<i>Replaces 05.003.137S</i>					
					No Tooth Offset

Compatible with DePuy® Synthes®  
Large Bone Oscillating Saw System



DEPUY® SYNTHES®

# BUSA® Large Bone Saw Blades - Oscillating Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System




DEPUY® SYNTHES®


Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM519-10</b>	mm	50	27	0.40	0.66
	in	1.97	1.06	0.016	0.026
<i>Replaces 519.10</i>					
					<i>Tooth Offset</i>
<b>KM519-15</b>	mm	50	14	0.40	0.66
	in	1.97	0.55	0.016	0.026
<i>Replaces 519.15</i>					
					<i>Tooth Offset</i>
<b>KM519-17</b>	mm	70	27	0.64	0.81
	in	2.76	1.06	0.025	0.032
<i>Replaces 519.17</i>					
					<i>Tooth Offset</i>
<b>KM519-21</b>	mm	70	27	0.99	1.19
	in	2.76	1.06	0.039	0.047
<i>Replaces 519.21</i>					
					<i>Tooth Offset</i>

# BUSA® Large Bone Saw Blades - Oscillating Replacements - DePuy® Synthes®

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM519-53SP</b>	mm	75	20	0.72	0.74
	in	2.95	0.79	0.028	0.029
<i>Replaces N/A</i>					
					<i>Tooth Offset</i>  <i>1 per Pack</i> <i>Non-Sterile</i>

## Special:

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-1201-84</b>	mm	86	19.5	1.27	1.27
	in	3.39	0.77	0.050	0.050
<i>Replaces N/A (Similar to Linvatec® 5071-201)</i>					
					<i>No Tooth Offset</i>  <i>Opposing Teeth</i> <i>Chromium Coating</i>

Compatible with DePuy® Synthes®  
Large Bone Oscillating Saw System



DEPUY® SYNTHES®