





BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-AM2U	Acorn	mm	7.5	9.4
		in	.30	.37
Replaces Linvatec® 7021-382 Replaces Medtronic® Midas Rex® AM-2 Replaces Stryker® 5130-30-75, 5230-30-75, 5430-30-75				
				
KM-M8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-237 Replaces Medtronic® Midas Rex® M-8 Replaces Stryker® 5120-107-30, 5220-107-30, 5420-107-30				
				
KM-AM8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-337 Replaces Medtronic® Midas Rex® AM-8 Replaces Stryker® 5130-107-30, 5230-107-30, 5430-107-30				
				
KM-AM8DU	Match Head Diamond	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-379 Replaces Medtronic® Midas Rex® AM-8D Replaces Stryker® 5130-107-130, 5230-107-130, 5430-107-32				
				

Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™







BUSA® High Speed Burs - Replacements

Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™





Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M51U	Metal Cutter (Diamond Wheel)	mm	25.4	0.8
		in	1.00	0.03
Replaces Linvatec® 7021-292 Replaces Medtronic® Midas Rex® M-51 Replaces Stryker® 5120-130-250				
				
KM-M52U	Metal Cutter (Carbide)	mm	3.1	17.5
		in	0.12	0.69
Replaces Linvatec® 7021-294 Replaces Medtronic® Midas Rex® M-52 Replaces Stryker® 5120-80-30, 5420-80-30				
				
KM-U295	Oval	mm	4.0	8.0
		in	0.16	0.31
Replaces Linvatec® 7021-295 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-15-140				
				
KM5110-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5110-15-40				
				

STRYKER®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM5130-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-15-40				
				
KM-M340DU	Round Diamond	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-264 Replaces Medtronic® Midas Rex® M-340D Replaces Stryker® 5120-12-40				
				
KM-U265	Round Diamond	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-265 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-12-50				
				
KM-M320U	Round Fluted	mm	2.0	2.0
		in	0.08	0.08
Replaces Linvatec® 7021-224 Replaces Medtronic® Midas Rex® M-320 Replaces Stryker® 5120-10-20, 5220-10-20				
				

Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™







STRYKER®
 Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Compatible with High Speed Systems: Linvatec® UltraPower®;
Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-226 Replaces Medtronic® Midas Rex® M-330 Replaces Stryker® 5120-10-30, 5220-10-30				
				
KM5120-10-230	Round Fluted Extended Length	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-230, 5220-10-230				
				
KM-M340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-227 Replaces Medtronic® Midas Rex® M-340 Replaces Stryker® 5120-10-40, 5220-10-40				
				
KM5120-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-240, 5220-10-240				
				

STRYKER®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-228 Replaces Medtronic® Midas Rex® M-35 Replaces Stryker® 5120-10-50, 5220-10-50, 5420-10-50				
				
KM5120-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-250, 5220-10-250				
				
KM-M33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-229 Replaces Medtronic® Midas Rex® M-33 Replaces Stryker® 5120-10-60, 5220-10-60				
				
KM-M32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-231 Replaces Medtronic® Midas Rex® M-32 Replaces Stryker® 5120-10-75, 5220-10-75				
				

Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™







STRYKER®
 Replacements

BUSA® High Speed Burs - Replacements

Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M31U	Round Fluted	mm	9.0	9.0
		in	0.35	0.35
Replaces Linvatec® 7021-232 Replaces Medtronic® Midas Rex® M-31 Replaces Stryker® 5120-10-90, 5220-10-90				
				
KM-AM330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-326 Replaces Medtronic® Midas Rex® AM-330 Replaces Stryker® 5130-10-30, 5230-10-30				
				
KM-AM340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-327 Replaces Medtronic® Midas Rex® AM-340 Replaces Stryker® 5130-10-40, 5230-10-40				
				
KM5130-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-240, 5230-10-240				
				

STRYKER®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-AM35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-328 Replaces Medtronic® Midas Rex® AM-35 Replaces Stryker® 5130-10-50, 5230-10-50, 5430-10-50				
				
KM5130-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-250, 5230-10-250				
				
KM-AM33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-329 Replaces Medtronic® Midas Rex® AM-33 Replaces Stryker® 5130-10-60, 5230-10-60				
				
KM-AM32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-331 Replaces Medtronic® Midas Rex® AM-32 Replaces Stryker® 5130-10-75, 5230-10-75				
				

Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



STRYKER®
 Replacements