





BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-AM2U	Acorn	mm	7.5	9.4
		in	.30	.37
Replaces Linvatec® 7021-382 Replaces Medtronic® Midas Rex® AM-2 Replaces Stryker® 5130-30-75, 5230-30-75, 5430-30-75				
				
KM-M8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-237 Replaces Medtronic® Midas Rex® M-8 Replaces Stryker® 5120-107-30, 5220-107-30, 5420-107-30				
				
KM-AM8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-337 Replaces Medtronic® Midas Rex® AM-8 Replaces Stryker® 5130-107-30, 5230-107-30, 5430-107-30				
				
KM-AM8DU	Match Head Diamond	mm	2.9	3.8
		in	0.11	0.15
Replaces Linvatec® 7021-379 Replaces Medtronic® Midas Rex® AM-8D Replaces Stryker® 5130-107-130, 5230-107-130, 5430-107-32				
				





Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™





Compatible with High Speed Systems: Linvatec® UltraPower®;
Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M51U	Metal Cutter (Diamond Wheel)	mm	25.4	0.8
		in	1.00	0.03
Replaces Linvatec® 7021-292 Replaces Medtronic® Midas Rex® M-51 Replaces Stryker® 5120-130-250				
				
KM-M52U	Metal Cutter (Carbide)	mm	3.1	17.5
		in	0.12	0.69
Replaces Linvatec® 7021-294 Replaces Medtronic® Midas Rex® M-52 Replaces Stryker® 5120-80-30, 5420-80-30				
				
KM-U295	Oval	mm	4.0	8.0
		in	0.16	0.31
Replaces Linvatec® 7021-295 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-15-140				
				
KM5110-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5110-15-40				
				

CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM5130-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-15-40				
				
KM-M340DU	Round Diamond	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-264 Replaces Medtronic® Midas Rex® M-340D Replaces Stryker® 5120-12-40				
				
KM-U265	Round Diamond	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-265 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-12-50				
				
KM-M320U	Round Fluted	mm	2.0	2.0
		in	0.08	0.08
Replaces Linvatec® 7021-224 Replaces Medtronic® Midas Rex® M-320 Replaces Stryker® 5120-10-20, 5220-10-20				
				





Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-226 Replaces Medtronic® Midas Rex® M-330 Replaces Stryker® 5120-10-30, 5220-10-30				
				
KM5120-10-230	Round Fluted Extended Length	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-230, 5220-10-230				
				
KM-M340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-227 Replaces Medtronic® Midas Rex® M-340 Replaces Stryker® 5120-10-40, 5220-10-40				
				
KM5120-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-240, 5220-10-240				
				



CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Compatible with High Speed Systems: Linvatec® UltraPower®, Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-228 Replaces Medtronic® Midas Rex® M-35 Replaces Stryker® 5120-10-50, 5220-10-50, 5420-10-50				
				
KM5120-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5120-10-250, 5220-10-250				
				
KM-M33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-229 Replaces Medtronic® Midas Rex® M-33 Replaces Stryker® 5120-10-60, 5220-10-60				
				
KM-M32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-231 Replaces Medtronic® Midas Rex® M-32 Replaces Stryker® 5120-10-75, 5220-10-75				
				







CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Compatible with High Speed Systems: Linvatec® UltraPower®;
Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™

Order Number	Description	Unit	Head Diameter	Head Length
KM-M31U	Round Fluted	mm	9.0	9.0
		in	0.35	0.35
Replaces Linvatec® 7021-232 Replaces Medtronic® Midas Rex® M-31 Replaces Stryker® 5120-10-90, 5220-10-90				
				
KM-AM330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Linvatec® 7021-326 Replaces Medtronic® Midas Rex® AM-330 Replaces Stryker® 5130-10-30, 5230-10-30				
				
KM-AM340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® 7021-327 Replaces Medtronic® Midas Rex® AM-340 Replaces Stryker® 5130-10-40, 5230-10-40				
				
KM5130-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-240, 5230-10-240				
				



CONMED® LINVATEC®
Replacements

BUSA® High Speed Burs - Replacements Linvatec®, Medtronic® Midas Rex® & Stryker®

Compatible with High Speed Systems: Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™


Order Number	Description	Unit	Head Diameter	Head Length
KM-AM35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® 7021-328 Replaces Medtronic® Midas Rex® AM-35 Replaces Stryker® 5130-10-50, 5230-10-50, 5430-10-50				
				
KM5130-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Linvatec® N/A Replaces Medtronic® Midas Rex® N/A Replaces Stryker® 5130-10-250, 5230-10-250				
				
KM-AM33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Linvatec® 7021-329 Replaces Medtronic® Midas Rex® AM-33 Replaces Stryker® 5130-10-60, 5230-10-60				
				
KM-AM32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Linvatec® 7021-331 Replaces Medtronic® Midas Rex® AM-32 Replaces Stryker® 5130-10-75, 5230-10-75				
				

Compatible with High Speed Systems: Linvatec® UltraPower®,
 Medtronic® Midas Rex® Classic I & II; Stryker® CORE & TPS™



BUSA® High Speed Burs - Replacements - Linvatec®

Compatible with Linvatec® UltraPower® High Speed Systems







Order Number	Description	Unit	Head Diameter	Head Length
KM-U815	Craniotome Bur	mm	2.2	16.0
		in	0.09	0.63
Replaces Linvatec® 7021-815 Replaces Medtronic® Midas Rex® N/A Replaces Stryker® N/A				
				

Compatible with Linvatec®
UltraPower® High Speed Systems

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C1803	Lindemann - Carbide	mm	51.0	1.9	13.0	4
		in	2.01	0.07	0.51	
Replaces N/A						
						
03-C0901	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces 091-136, 5091-236						
						
03-C0905	Oval - Carbide	mm	48.0	5.5	10.2	10
		in	1.89	0.22	0.40	
Replaces 5091-138, 5091-238						
						
09-C0901	Oval - Carbide	mm	70.0	4.0	8.0	8
		in	2.76	0.16	0.31	
Replaces 5092-136, 5092-236						
						
09-C0905	Oval - Carbide	mm	70.0	5.5	10.2	10
		in	2.76	0.22	0.40	
Replaces 5092-138, 5092-238						
						
12-C0901	Oval - Carbide	mm	95.0	4.0	8.0	8
		in	3.74	0.16	0.31	
Replaces 5093-136, 5093-236						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements







BUSA® Micro Burs - Standard & Oral Surgery

2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-C0905	Oval - Carbide	mm	95.0	5.5	10.2	10
		in	3.74	0.22	0.40	
Replaces 5093-138, 5093-238						
						
04-C0101	Round - Carbide	mm	51.0	1	1	6
		in	2.01	0.04	0.04	
Replaces 5091-120, 5091-220						
						
04-C0103	Round - Carbide	mm	51.0	1.5	1.5	6
		in	2.01	0.06	0.06	
Replaces 5091-122, 5091-222						
						
04-C0104	Round - Carbide	mm	51.0	1.8	1.8	6
		in	2.01	0.07	0.07	
Replaces 5091-121, 5091-221						
						
04-C0105	Round - Carbide	mm	51.0	2.0	2.0	6
		in	2.01	0.08	0.08	
Replaces 5091-124, 5091-224						
						
04-C0106	Round - Carbide	mm	51.0	2.3	2.3	8
		in	2.01	0.09	0.09	
Replaces 5091-123, 5091-125, 5091-223, 5091-225						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0107	Round - Carbide	mm	51.0	3.0	3.0	8
		in	2.01	0.12	0.12	
Replaces 5091-126, 5091-226						
						
04-C0108	Round - Carbide	mm	51.0	4.0	4.0	8
		in	2.01	0.16	0.16	
Replaces 5091-128, 5091-228						
						
04-C0109	Round - Carbide	mm	51.0	5.0	5.0	10
		in	2.01	0.20	0.20	
Replaces 5091-130, 5091-230						
						
04-C0110	Round - Carbide	mm	51.0	6.0	6.0	8
		in	2.01	0.24	0.24	
Replaces M60R						
						
04-C0170	Round - Carbide	mm	51.0	6.5	6.5	10
		in	2.01	0.26	0.26	
Replaces 5091-132, 5091-232						
						
09-C0109	Round - Carbide	mm	70.0	5.0	5.0	10
		in	2.76	0.20	0.20	
Replaces 5092-130, 5092-230						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-C0110	Round - Carbide	mm	70.0	6.0	6.0	8
		in	2.76	0.24	0.24	
Replaces L60R						
						
09-C0170	Round - Carbide	mm	70.0	6.5	6.5	10
		in	2.76	0.26	0.26	
Replaces 5092-132, 5092-232						
						
10-C0101	Round - Carbide	mm	71.5	1.0	1.0	6
		in	2.81	0.04	0.04	
Replaces 5092-120, 5092-220						
						
10-C0103	Round - Carbide	mm	72.0	1.5	1.5	6
		in	2.83	0.06	0.06	
Replaces 5092-122, 5092-222						
						
10-C0105	Round - Carbide	mm	72.5	2.0	2.0	6
		in	2.85	0.08	0.08	
Replaces 5092-124, 5092-224						
						
10-C0107	Round - Carbide	mm	73.0	3.0	3.0	8
		in	2.87	0.12	0.12	
Replaces 5092-126, 5092-226						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0108	Round - Carbide	mm	73.0	4.0	4.0	8
		in	2.87	0.16	0.16	
Replaces 5092-128, 5092-228						
						
12-C0105	Round - Carbide	mm	95.0	2.0	2.0	6
		in	3.74	0.08	0.08	
Replaces 5093-124, 5093-224						
						
12-C0108	Round - Carbide	mm	95.0	4.0	4.0	8
		in	3.74	0.16	0.16	
Replaces 5093-128, 5093-228						
						
12-C0109	Round - Carbide	mm	95.0	5.0	5.0	10
		in	3.74	0.20	0.20	
Replaces 5093-130, 5093-230						
						
12-C0170	Round - Carbide	mm	95.0	6.5	6.5	10
		in	3.74	0.26	0.26	
Replaces 5093-132, 5093-232						
						
04-C0601	Tapered - Carbide	mm	51.0	0.9	5.2	6
		in	2.01	0.035	0.20	
Replaces 5091-115, 5091-215						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0604	Tapered - Carbide	mm	51.0	1.0	6.0	6
		in	2.01	0.04	0.24	
Replaces N/A						
						
04-C0606	Tapered - Carbide	mm	51.0	1.2	4.2	6
		in	2.01	0.05	0.17	
Replaces N/A						
						
04-C0623	Tapered - Carbide	mm	51.0	1.2	6.0	6
		in	2.01	0.05	0.24	
Replaces N/A						
						
04-C0609	Tapered - Carbide	mm	51.0	1.5	10.0	6
		in	2.01	0.06	0.39	
Replaces 5091-102, 5091-202						
						
04-C0607	Tapered - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces 5091-118, 5091-218						
						
04-C0617	Tapered - Carbide	mm	51.0	2.0	7.5	6
		in	2.01	0.08	0.30	
Replaces 5091-106, 5091-206						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0615	Tapered - Carbide	mm	51.0	2.1	4.8	6
		in	2.01	0.08	0.19	
Replaces N/A						
						
04-C0612	Tapered - Carbide	mm	51.0	2.1	4.9	6
		in	2.01	0.08	0.19	
Replaces 5091-113, 5091-213						
						
06-C0619	Tapered - Carbide	mm	57.0	2.0	19.0	6
		in	2.24	0.08	0.75	
Replaces 5091-114, 5091-214						
						
09-C0613	Tapered - Carbide	mm	70.0	1.73	12.0	6
		in	2.76	0.07	0.47	
Replaces 5092-103, 5092-203						
						
10-C0617	Tapered - Carbide	mm	73.0	2.0	7.5	6
		in	2.87	0.08	0.30	
Replaces 5092-106, 5092-206						
						
11-C0619	Tapered - Carbide	mm	80.0	2.0	19.0	6
		in	3.15	0.08	0.75	
Replaces 5092-114, 5092-214						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
03-C0312	Straight - Cross-Cut - Carbide	mm	48.0	4.0	10.0	6
		in	1.89	0.16	0.39	
Replaces 5091-154, 5091-254						
						
04-C0302	Straight - Cross-Cut - Carbide	mm	51.0	1.0	5.0	6
		in	2.01	0.04	0.20	
Replaces 5091-101, 5091-201						
						
04-C0404	Straight - Cross-Cut - Carbide	mm	51.0	1.2	3.7	6
		in	2.01	0.05	0.15	
Replaces 5091-119, 5091-219						
						
04-C0306	Straight - Cross-Cut - Carbide	mm	51.0	1.5	7.5	6
		in	2.01	0.06	0.30	
Replaces 5091-104, 5091-204						
						
04-C0307	Straight - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces 5091-117, 5091-217						
						
04-C0310	Straight - Cross-Cut - Carbide	mm	51.0	2.04	10.0	8
		in	2.01	0.08	0.39	
Replaces 5091-152, 5091-252						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0306	Straight - Carbide	mm	73.0	1.5	7.5	6
		in	2.87	0.06	0.30	
Replaces 5092-104, 5092-204						
						
10-C0310	Straight - Cross-Cut - Carbide	mm	73.0	2.04	10.0	8
		in	2.87	0.08	0.39	
Replaces 5092-152, 5092-252						
						
03-D0911*	Oval - Diamond - Dermabrader	mm	48.0	4.0	8.0	Medium
		in	1.89	0.16	0.31	
Replaces 5091-170, 5091-270						
						
*Check Manufacturer's recommended bur guards and speed for handpiece being used.						
03-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
Replaces 5091-164, 5091-264						
						
03-D0138	Round - Diamond	mm	48.0	6.0	6.0	Coarse
		in	1.89	0.24	0.24	
Replaces 5091-168, 5091-268						
						
04-D0139	Round - Diamond	mm	51.0	2.0	2.0	Medium
		in	2.01	0.08	0.08	
Replaces 5091-160, 5091-260						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-D0133	Round - Diamond	mm	51.0	3.0	3.0	Medium
		in	2.01	0.12	0.12	
Replaces 5091-162, 5091-262						
						
09-D0128	Round - Diamond	mm	70.0	1.0	1.0	Fine
		in	2.76	0.04	0.04	
Replaces N/A						
						
09-D0140	Round - Diamond	mm	70.0	2.3	2.3	Medium
		in	2.76	0.09	0.09	
Replaces N/A						
						
09-D0132	Round - Diamond	mm	70.0	2.9	2.9	Medium
		in	2.76	0.11	0.11	
Replaces N/A						
						
09-D0135	Round - Diamond	mm	70.0	4.0	4.0	Medium
		in	2.76	0.16	0.16	
Replaces 5092-164, 5092-264						
						
09-D0141	Round - Diamond	mm	70.0	5.0	5.0	Coarse
		in	2.76	0.20	0.20	
Replaces 5092-166, 5092-266						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-D0138	Round - Diamond	mm	70.0	6.0	6.0	Coarse
		in	2.76	0.24	0.24	
Replaces 5092-168, 5092-268						
						
10-D0139	Round - Diamond	mm	73.0	2.0	2.0	Medium
		in	2.87	0.08	0.08	
Replaces 5092-160, 5092-260						
						
10-D0133	Round - Diamond	mm	73.0	3.0	3.0	Medium
		in	2.87	0.12	0.12	
Replaces 5092-162, 5092-262						
						
12-D0133	Round - Diamond	mm	95.0	3.0	3.0	Medium
		in	3.74	0.12	0.12	
Replaces 5093-162, 5093-262						
						
12-D0135	Round - Diamond	mm	95.0	4.0	4.0	Medium
		in	3.74	0.16	0.16	
Replaces 5093-164, 5093-264						
						
12-D0141	Round - Diamond	mm	95.0	5.0	5.0	Coarse
		in	3.74	0.20	0.20	
Replaces 5093-166, 5093-266						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
08-S1802	Lindemann - Stainless Steel	mm	67.0	2.1	25.0	4
		in	2.64	0.083	.98	
Replaces 5091-112, 5091-212						
						
03-S1101	Flame - Stainless Steel	mm	48.0	6.0	13.0	14
		in	1.89	0.24	0.51	
Replaces N/A						
						
09-S0101	Round - Stainless Steel	mm	70.0	1.0	1.0	6
		in	2.76	0.04	0.04	
Replaces N/A						
						
09-S0150	Round - Stainless Steel	mm	70.0	1.4	1.4	8
		in	2.76	0.058	0.058	
Replaces N/A						
						
09-S0151	Round - Stainless Steel	mm	70.0	1.8	1.8	8
		in	2.76	0.07	0.07	
Replaces N/A						
						
09-S0115	Round - Stainless Steel	mm	70.0	2.4	2.4	8
		in	2.76	0.091	0.091	
Replaces N/A						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-S0153	Round - Stainless Steel	mm	70.0	3.1	3.1	10
		in	2.76	0.12	0.12	
<i>Replaces N/A</i>						
						
09-S0120	Round - Stainless Steel	mm	70.0	4.0	4.0	12
		in	2.76	0.16	0.16	
<i>Replaces N/A</i>						
						
09-S0155	Round - Stainless Steel	mm	70.0	5.0	5.0	12
		in	2.76	0.20	0.20	
<i>Replaces N/A</i>						
						
09-S0156	Round - Stainless Steel	mm	70.0	6.0	6.0	14
		in	2.76	0.24	0.24	
<i>Replaces N/A</i>						
						
09-S0125	Round - Stainless Steel	mm	70.0	7.0	7.0	16
		in	2.76	0.28	0.28	
<i>Replaces N/A</i>						
						
12-S0153	Round - Stainless Steel	mm	95.0	3.1	3.1	10
		in	3.74	0.12	0.12	
<i>Replaces N/A</i>						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems







Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®, MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-S0123	Round - Stainless Steel	mm	95.0	4.5	4.5	12
		in	3.74	0.179	0.18	
Replaces N/A						
						
12-S0155	Round - Stainless Steel	mm	95.0	5.0	5.0	12
		in	3.74	0.20	0.20	
Replaces N/A						
						
04-S1501	Shannon Style - Tapered - Stainless Steel	mm	51.0	2.3	12.0	4
		in	2.01	0.09	0.47	
Replaces N/A						
						
04-S1502	Shannon Style - Straight - Stainless Steel	mm	51.0	2.4	13.2	3
		in	2.01	0.091	0.52	
Replaces N/A						
						
06-S1601	Swanson Reamer - Stainless Steel	mm	57.0	2.0	10.0	5
		in	2.24	0.08	0.39	
Replaces 5091-142, 5091-242						
						
07-S1602	Swanson Reamer - Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces 5091-143, 5091-243						
						

BUSA® Micro Burs - Standard & Oral Surgery 2.35 mm Straight Shank - Replacements - Linvatec®

Compatible with Micro Drills: Linvatec® MicroPower®, MicroChoice®, Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills; All BUSA® Small Bone Drill Systems





Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
07-S1605	Swanson Reamer - Stainless Steel	mm	62.0	4.0	16.5	6
		in	2.44	0.16	0.65	
Replaces 5091-144, 5091-244						
						
01-S1703	Wire Pass Drill - Stainless Steel	mm	37.0	1.5	19.0	2
		in	1.46	0.059	0.75	
Replaces 5090-148, 5090-248						
						
01-S1704	Wire Pass Drill - Stainless Steel	mm	37.0	2.0	19.0	2
		in	1.46	0.079	0.75	
Replaces 5090-149, 5090-249						
						
04-S1701	Wire Pass Drill - Stainless Steel	mm	51.0	1.0	19.0	2
		in	2.01	0.039	0.75	
Replaces 5091-147, 5091-247						
						
04-S1703	Wire Pass Drill - Stainless Steel	mm	51.0	1.5	19.0	2
		in	2.01	0.059	0.75	
Replaces 5091-148, 5091-248						
						
04-S1704	Wire Pass Drill - Stainless Steel	mm	51.0	2.0	19.0	2
		in	2.01	0.079	0.75	
Replaces 5091-149, 5091-249						
						

Compatible with Micro Drills: Linvatec® - MicroPower®, MicroChoice®,
Micro 100™, Micro E® & Surgairtome Two®; MicroAire® Small Bone Drills;
All BUSA® Small Bone Drill Systems

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Cebotome® Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Cebotome®







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
13-S01189	Round - Stainless Steel	mm	30.1	5.0	5.0	6
		in	1.19	0.20	0.20	
<i>Replaces 5052-189</i>						
						
13-S01184	Round - Stainless Steel	mm	30.3	7.0	7.0	10
		in	1.19	0.28	0.28	
<i>Replaces 5052-184</i>						
						
13-S01185	Round - Stainless Steel	mm	31.5	8.0	8.0	10
		in	1.24	0.31	0.31	
<i>Replaces 5052-185</i>						
						
13-S01186	Round - Stainless Steel	mm	32.6	9.0	9.0	10
		in	1.28	0.35	0.35	
<i>Replaces 5052-186</i>						
						

Compatible with Linvatec®
Drill: Cebotome®

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
70-C0108	Round - Carbide	mm	30.0	4.0	4.0	8
		in	1.18	0.16	0.16	
<i>Replaces 5056-19, 5056-119</i>						
						
70-C0109	Round - Carbide	mm	30.0	5.0	5.0	10
		in	1.18	0.20	0.20	
<i>Replaces 5056-18, 5056-118</i>						
						
70-C0170	Round - Carbide	mm	30.0	6.5	6.5	10
		in	1.18	0.26	0.26	
<i>Replaces 5056-17, 5056-117</i>						
						
71-C0107	Round - Carbide	mm	38.0	3.0	3.0	8
		in	1.50	0.12	0.12	
<i>Replaces 5056-28, 5056-128</i>						
						
71-C0108	Round - Carbide	mm	38.0	4.0	4.0	8
		in	1.50	0.16	0.16	
<i>Replaces 5056-27, 5056-127</i>						
						
71-C0109	Round - Carbide	mm	38.0	5.0	5.0	10
		in	1.50	0.20	0.20	
<i>Replaces 5056-26, 5056-126</i>						
						

Compatible with Linvatec®
Drill: Osteon & Ototome


CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
71-C0170	Round - Carbide	mm	38.0	6.5	6.5	10
		in	1.50	0.26	0.26	
Replaces 5056-25, 5056-125						
						
72-C0107	Round - Carbide	mm	48.0	3.0	3.0	8
		in	1.89	0.12	0.12	
Replaces 5056-46, 5056-146						
						
72-C0108	Round - Carbide	mm	48.0	4.0	4.0	8
		in	1.89	0.16	0.16	
Replaces 5056-45, 5056-145						
						
72-C0109	Round - Carbide	mm	48.0	5.0	5.0	10
		in	1.89	0.20	0.20	
Replaces 5056-44, 5056-144						
						
75-S0101	Round - Stainless Steel	mm	48.0	1.0	1.0	6
		in	1.89	0.04	0.04	
Replaces 5056-64, 5056-164						
						
75-S0148	Round - Stainless Steel	mm	48.0	0.56	0.5	6
		in	1.89	0.022	0.02	
Replaces 5056-66, 5056-166						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome







Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
75-S0149	Round - Stainless Steel	mm	48.0	0.74	0.74	6
		in	1.89	0.03	0.03	
Replaces 5056-65, 5056-165						
						
75-S0150	Round - Stainless Steel	mm	48.0	1.4	1.4	8
		in	1.89	0.058	0.058	
Replaces 5056-40, 5056-140						
						
75-S0151	Round - Stainless Steel	mm	48.0	1.8	1.8	8
		in	1.89	0.07	0.07	
Replaces 5056-42, 5056-142						
						
71-D0135	Round - Diamond	mm	38.0	4.0	4.0	Medium
		in	1.50	0.16	0.16	
Replaces 5056-48, 5056-148						
						
71-D0141	Round - Diamond	mm	38.0	5.0	5.0	Coarse
		in	1.50	0.20	0.20	
Replaces 5056-50, 5056-150						
						
72-D0133	Round - Diamond	mm	48.0	3.0	3.0	Medium
		in	1.89	0.12	0.12	
Replaces 5056-59, 5056-159						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements

BUSA® Micro Burs - Standard Osteon/Ototome Shank - Replacements - Linvatec®

Compatible with Linvatec® Drill: Osteon & Ototome

Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
72-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
<i>Replaces 5056-52, 5056-152</i>						
						
72-D0141	Round - Diamond	mm	48.0	5.0	5.0	Coarse
		in	1.89	0.20	0.20	
<i>Replaces N/A</i>						
						
75-D0128	Round - Diamond	mm	48.0	1.0	1.0	Fine
		in	1.89	0.04	0.04	
<i>Replaces 5056-58, 5056-158</i>						
						
75-D0139	Round - Diamond	mm	48.0	2.0	2.0	Medium
		in	1.89	0.08	0.08	
<i>Replaces 5056-62, 5056-162</i>						
						
75-D0142	Round - Diamond	mm	48.0	0.54	0.54	X-Fine
		in	1.89	0.021	0.021	
<i>Replaces 5056-61, 5056-161</i>						
						
75-D0143	Round - Diamond	mm	48.0	0.72	0.72	Fine
		in	1.89	0.03	0.03	
<i>Replaces 5056-60, 5056-160</i>						
						

Compatible with Linvatec®
Drill: Osteon & Ototome

CONMED® LINVATEC®
Replacements