






BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thin) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-51	mm	18.5	6.5	0.38	0.43
	in	0.73	0.26	0.015	0.017
Replaces 2296-3-51					
 Thin					
KM3-101	mm	10.0	9.0	0.38	0.43
	in	0.39	0.35	0.015	0.017
Replaces 2296-3-101					
 Thin					
KM3-102	mm	18.5	16.5	0.38	0.43
	in	0.73	0.65	0.015	0.017
Replaces 2296-3-102					
 Thin					
KM3-103	mm	15.0	7.0	0.38	0.43
	in	0.59	0.28	0.015	0.017
Replaces 2296-3-103					
 Thin					
KM3-104	mm	34.5	16.5	0.38	0.43
	in	1.36	0.65	0.015	0.017
Replaces 2296-3-104					
 Thin					

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2



STRYKER®
Replacements






BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thin) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Compatible with Stryker® Sagittal & Oscillating
Saw Systems: RemB™, CORE, TPS™ & Command 2








STRYKER®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-105	mm	18.5	9.0	0.38	0.43
	in	0.73	0.35	0.015	0.017
Replaces 2296-3-105					
 Thin					
KM3-106	mm	34.5	13.0	0.38	0.43
	in	1.36	0.51	0.015	0.017
Replaces 2296-3-106					
 Thin					
KM3-108	mm	39.0	13.0	0.38	0.43
	in	1.54	0.51	0.015	0.017
Replaces 2296-3-108					
 Thin					
KM3-111	mm	25.0	9.0	0.38	0.43
	in	0.98	0.35	0.015	0.017
Replaces 2296-3-111					
 Thin					
KM3-114	mm	18.5	7.0	0.38	0.43
	in	0.73	0.28	0.015	0.017
Replaces 2296-3-114					
 Thin					

BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thin) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-115	mm	29.5	7.0	0.38	0.43
	in	1.16	0.28	0.015	0.017
Replaces 2296-3-115					
 Thin					
KM3-125	mm	31.0	9.0	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces 2296-3-125					
 Thin					
KM3-206	mm	34.5	12.0	0.38	0.43
	in	1.36	0.47	0.015	0.017
Replaces 2296-3-206					
 Thin					
KM3-212	mm	18.0	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces 2296-3-212					
 Thin					
KM3-225	mm	31.0	9.0	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces 2296-3-225					
 Thin					

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2



STRYKER®
Replacements






BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thin) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Compatible with Stryker® Sagittal & Oscillating
Saw Systems: RemB™, CORE, TPS™ & Command 2






STRYKER®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-254	mm	41.0	19.5	0.38	0.68
	in	1.61	0.77	0.015	0.027
Replaces 2296-3-254					
					Thin
KM3-410	mm	11.5	5.5	0.38	0.43
	in	0.45	0.22	0.015	0.017
Replaces 2296-3-410					
					Thin
KM3-411	mm	11.5	9.0	0.38	0.43
	in	0.45	0.35	0.015	0.017
Replaces 2296-3-411					
					Thin
KM3-412	mm	18.0	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces 2296-3-412					
					Thin
KM3-413	mm	18.0	2.5	0.38	0.43
	in	0.71	0.10	0.015	0.017
Replaces 2296-3-413					
					Thin

BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thin) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM3-414	mm	25.0	5.5	0.38	0.43
	in	0.98	0.22	0.015	0.017
Replaces 2296-3-414					
					Thin
KM3-511	mm	25.0	9.0	0.38	0.64
	in	0.98	0.35	0.015	0.025
Replaces 2296-3-511					
					Aggressive Cut Thin
KM3-525	mm	31.0	9.0	0.38	0.64
	in	1.22	0.35	0.015	0.025
Replaces 2296-3-525					
					Aggressive Cut Thin

Compatible with Stryker® Sagittal & Oscillating
Saw Systems: RemB™, CORE, TPS™ & Command 2



STRYKER®
Replacements






BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thick) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Compatible with Stryker® Sagittal & Oscillating
Saw Systems: RemB™, CORE, TPS™ & Command 2








STRYKER®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-101	mm	10.0	9.0	0.51	0.56
	in	0.39	0.35	0.020	0.022
Replaces 2296-33-101					
 <p style="text-align: right;">Thick</p>					
KM33-102	mm	18.5	16.5	0.51	0.56
	in	0.73	0.65	0.020	0.022
Replaces 2296-33-102					
 <p style="text-align: right;">Thick</p>					
KM33-103	mm	15.0	7.00	0.51	0.56
	in	0.59	0.28	0.020	0.022
Replaces 2296-33-103					
 <p style="text-align: right;">Thick</p>					
KM33-104	mm	34.5	16.5	0.51	0.56
	in	1.36	0.65	0.020	0.022
Replaces 2296-33-104					
 <p style="text-align: right;">Thick</p>					
KM33-105	mm	18.5	9.0	0.51	0.56
	in	0.73	0.35	0.020	0.022
Replaces 2296-33-105					
 <p style="text-align: right;">Thick</p>					

BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thick) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-106	mm	34.5	13.0	0.51	0.56
	in	1.36	0.51	0.020	0.022
Replaces 2296-33-106					
 <p style="text-align: right;">Thick</p>					
KM33-108	mm	39.0	13.0	0.51	0.56
	in	1.54	0.51	0.020	0.022
Replaces 2296-33-108					
 <p style="text-align: right;">Thick</p>					
KM33-111	mm	25.0	9.0	0.51	0.56
	in	0.98	0.35	0.020	0.022
Replaces 2296-33-111					
 <p style="text-align: right;">Thick</p>					
KM33-114	mm	18.5	7.0	0.51	0.56
	in	0.73	0.28	0.020	0.022
Replaces 2296-33-114					
 <p style="text-align: right;">Thick</p>					
KM33-115	mm	29.5	7.0	0.51	0.56
	in	1.16	0.28	0.020	0.022
Replaces 2296-33-115					
 <p style="text-align: right;">Thick</p>					

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2








BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thick) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Compatible with Stryker® Sagittal & Oscillating
Saw Systems: RemB™, CORE, TPS™ & Command 2








STRYKER®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-125	mm	31.0	9.0	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces 2296-33-125					
 Thick					
KM33-206	mm	34.5	12.0	0.51	0.56
	in	1.36	0.47	0.020	0.022
Replaces 2296-33-206					
 Thick					
KM33-212	mm	18.0	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces 2296-33-212					
 Thick					
KM33-225	mm	31.0	9.0	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces 2296-33-225					
 Thick					
KM33-410	mm	11.5	5.5	0.51	0.56
	in	0.45	0.22	0.020	0.022
Replaces 2296-33-410					
 Thick					

BUSA® Small Bone Sagittal & Oscillating Saw Blades (Thick) Replacements - Stryker®

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM33-411	mm	11.5	9.0	0.51	0.56
	in	0.45	0.35	0.020	0.022
Replaces 2296-33-411					
 Thick					
KM33-412	mm	18.0	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces 2296-33-412					
 Thick					
KM33-413	mm	18.0	2.5	0.51	0.56
	in	0.71	0.10	0.020	0.022
Replaces 2296-33-413					
 Thick					
KM33-414	mm	25.0	5.5	0.51	0.56
	in	0.98	0.22	0.020	0.022
Replaces 2296-33-414					
 Thick					
KM33-521	mm	25.0	9.0	0.51	0.74
	in	0.98	0.35	0.020	0.029
Replaces 2296-33-521					
 Thick					

Compatible with Stryker® Sagittal & Oscillating Saw Systems: RemB™, CORE, TPS™ & Command 2







STRYKER®
Replacements

BUSA® Small Bone Reciprocating Saw Blades Replacements - Stryker®

Compatible with Stryker® Reciprocating Saw Systems: RemB™, CORE, TPS™ & Command 2





Compatible with Stryker® Reciprocating Saw
Systems: RemB™, CORE, TPS™ & Command 2

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM5100-137-107	mm	23.0	6.6	0.38	0.51
	in	0.91	0.26	0.015	0.020
<i>Replaces 5100-137-107</i>					
					<i>Thin</i>
KM5100-137-123	mm	27.0	6.4	0.38	0.51
	in	1.06	0.25	0.015	0.020
<i>Replaces 5100-137-123</i>					
					<i>Thin</i>
KM5100-137-133	mm	21.8	6.4	0.38	0.51
	in	0.86	0.25	0.015	0.020
<i>Replaces 5100-137-133</i>					
					<i>Thin</i>
KM5100-137-233	mm	27.4	6.0	0.38	0.51
	in	1.08	0.24	0.015	0.020
<i>Replaces 5100-137-233</i>					
					<i>Thin</i>

STRYKER®
Replacements

BUSA® Small Bone Reciprocating Rasps - Replacements - Stryker®

Compatible with Stryker® Reciprocating Saw Systems: RemB™, CORE, TPS™ & Command 2



Order Number	Unit	Length	Max. Cut Width
KM5100-37-113	mm	11.0	5.5
	in	0.43	0.22
<i>Replaces 5100-37-113</i>			
			<i>Small Tear Cross-Cut</i>
KM5100-37-114	mm	13.0	6.5
	in	0.51	0.26
<i>Replaces 5100-37-114</i>			
			<i>Large Tear Cross-Cut</i>
KM5100-37-115	mm	11.0	5.5
	in	0.43	0.22
<i>Replaces 5100-37-115</i>			
			<i>Small Tear</i>
KM5100-37-116	mm	13.0	6.5
	in	0.51	0.26
<i>Replaces 5100-37-116</i>			
			<i>Large Tear</i>

Compatible with Stryker® Reciprocating Saw Systems: RemB™, CORE, TPS™ & Command 2



BUSA® Small Bone Oscillating Saw Blades Intra-Oral - Replacements - Stryker®

Compatible with Stryker® Oscillating Saw System: TPS™

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-2232	mm	7.6	11.7	0.38	0.61
	in	0.30	0.46	0.015	0.024
<i>Replaces 2296-31-32</i>					
 <p style="text-align: right;"><i>Thin</i></p>					
KM-2234	mm	12.0	11.5	0.38	0.61
	in	0.47	0.45	0.015	0.024
<i>Replaces 2296-31-34</i>					
 <p style="text-align: right;"><i>Thin</i></p>					






Compatible with Stryker®
Oscillating Saw System: TPS™



STRYKER®
Replacements

BUSA® Small Bone Oscillating Saw Blades Replacements - Stryker®

Compatible with Stryker® Oscillating Saw System: Micro



Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-103S	mm	25.5	9.5	0.38	0.64
	in	1.00	0.37	0.015	0.025
Replaces 277-31-103					
					
KM-410S	mm	12.5	5.5	0.38	0.64
	in	0.48	0.22	0.015	0.025
Replaces 277-31-410					
					
KM-412S	mm	20.5	5.5	0.38	0.64
	in	0.80	0.22	0.015	0.025
Replaces 277-31-412					
					
KM-413S	mm	18.5	9.5	0.38	0.64
	in	0.73	0.37	0.015	0.025
Replaces 277-31-413					
					
KM-414S	mm	25.5	5.5	0.38	0.64
	in	1.00	0.22	0.015	0.025
Replaces 277-31-414					
					

Compatible with Stryker®
Oscillating Saw System: Micro



BUSA® Small Bone Reciprocating Saw Blades Replacements - Stryker®

Compatible with Stryker® Reciprocating Saw System: Micro




Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-062	mm	27.2	6.0	0.64	1.00
	in	1.07	0.24	0.025	0.039
<i>Replaces 1675-133</i>					
					<i>Thick</i>
KM-062PT	mm	21.8	6.4	0.38	0.51
	in	0.86	0.25	0.015	0.020
<i>Replaces 296-37-133</i>					
					<i>Thin</i>

Compatible with Stryker®
Reciprocating Saw System: Micro

STRYKER®
Replacements

BUSA® Small Bone Reciprocating Rasps - Replacements - Stryker®

Compatible with Stryker® Reciprocating Saw System: Micro

Order Number	Unit	Length	Max. Cut Width
KM-800	mm	13.0	6.5
	in	0.51	0.26
<i>Replaces 1675-116</i>			
 <p style="text-align: right;"><i>Large Tear</i></p>			
KM-801	mm	11.0	5.5
	in	0.43	0.22
<i>Replaces 1675-115</i>			
 <p style="text-align: right;"><i>Small Tear</i></p>			
KM-802	mm	11.0	5.5
	in	0.43	0.22
<i>Replaces 1675-113</i>			
 <p style="text-align: right;"><i>Small Tear Cross-Cut</i></p>			

Compatible with Stryker®
Reciprocating Saw System: Micro

