







BUSA® Sternum Saw Blades - Replacements - Stryker®

Compatible with Stryker® Sternum Saw Systems: 7, 6, 5, 4, 2000, EHD & OrthoPower


Compatible with Stryker® Sternum Saw Systems:
7, 6, 5, 4, 2000, EHD & OrthoPower

STRYKER®
Replacements

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
ST6207-97-101	mm	30.5	6.50	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces 6207-97-101					
					Sterile
ST6207-97-101N	mm	30.5	6.50	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces 6207-97-101					
					Non-Sterile
KM-97	mm	32.0	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces 298-97-100					
					Sterile
KM-97N	mm	32.0	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces 298-97-100					
					Non-Sterile
KM-98	mm	32.0	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces 298-97-101					
					Sterile
KM-98N	mm	32.0	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces 298-97-101					
					Non-Sterile

BUSA® Sternum Saw Blades - Replacements - Stryker®

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR2108-125	mm	45.0	42.0	0.64	1.02
	in	1.77	1.64	0.025	0.040
Replaces 2108-125					
					

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD




BUSA® Sternum Saw Blades - Replacements - Stryker®

Compatible with Stryker® Sagittal Saw System: OrthoPower

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-325	mm	45.0	42.0	0.64	1.02
	in	1.77	1.64	0.025	0.040

Replaces 277-88-125



Sterile

Compatible with Stryker® Sagittal
Saw System: OrthoPower



STRYKER®
Replacements